

Station 2. "Therapeutic dentistry. Practical skills 2. »

Algorithm of action of the student on practical skill: "Total etching of hard tissues of 36 teeth with 37% gel of phosphoric acid".

1. Say hello.
2. Wash hands thoroughly and treat with antiseptic.
3. Wear a mask and glasses.
4. Wear gloves.
5. Explain the purpose of dental manipulation (orally).
6. Take a tube with 37% phosphoric acid in your right hand.
7. Replace the syringe cap with a cannula.
8. Dry the surface of 36 tooth.
9. Apply etchant on the prepared surface for 15-30 seconds.
10. Wash the etchant continuously for 15-30 seconds.
11. Do the control of etching (orally give explanation).
12. Remove gloves and place it to a container for used materials.
13. Wash hands thoroughly.

